P-9801

<https://www.gigahertz-optik.com/de-de/produkt/p-9801/>

Produkt-Tags: Dosimeter , Multikanal

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Messverfahren	<table><tr><th>Menüpunkt</th><th>Untermenüpunkt</th><th>Funktion</th></tr><tr><td rowspan="5">Display Mode</td><td>Absolute</td><td>displays the measurement respective of any offset and calibration factors programmed</td></tr><tr><td>Relative Percent</td><td>displays the measurement as percentage of a reference value</td></tr><tr><td>Relative Factor</td><td>displays the measurement referred to a reference value</td></tr><tr><td>Relative Logarithmic Dose</td><td>displays the measurement in dB in relation to a reference value</td></tr><tr><td>No Display</td><td>accumulates the single readings and displays the result as exposure for measured quantity</td></tr><tr><td>Range</td><td></td><td>no display of the selected channel</td></tr><tr><td rowspan="4">Detector</td><td>Ampere</td><td>selects the range 1 to 8 or switches into the autorange mode (9)</td></tr><tr><td>Manual</td><td>shows the reading in ampere units</td></tr><tr><td>Select Table Entry</td><td>allows entering a calibration factor manually</td></tr><tr><td>Select Detector Recalibration</td><td>selects calibration data table entry</td></tr><tr><td>Offset</td><td></td><td>selects calibration data table entry</td></tr><tr><td>Reference</td><td></td><td>performs recalibration of selected calibration table entry with a reference lamp</td></tr><tr><td rowspan="5">Main Mode</td><td>CW</td><td>performs an automatic offset adjustment</td></tr><tr><td>Color</td><td>sets a reference value, used in different display modes</td></tr><tr><td>Pulse</td><td>measurement display of single channels</td></tr><tr><td>Pulse Color</td><td>measurement display color mode</td></tr><tr><td>RS232</td><td>measurement of short and single light pulses</td></tr><tr><td rowspan="8">Setup</td><td>IEEE488</td><td>measurement of short and single light pulses, result displayed as color values</td></tr><tr><td>Measurement Time</td><td>remote control RS232</td></tr><tr><td>Analog Out</td><td>remote control IEEE488</td></tr><tr><td>Dose Run Time</td><td>sets the integration time</td></tr><tr><td>Dose Maximum</td><td>selects the channel switched to analog output</td></tr><tr><td>Display Illumination</td><td>sets the maximum time duration for exposure measurement</td></tr><tr><td>Display Digits</td><td>sets the max. dose level for the exposure measurement mode</td></tr><tr><td>Color Setup</td><td>display backlight on / off</td></tr><tr><td></td><td>IEEE488 Address</td><td>adjusts the number of displayed digits (4, 5, or automatic)</td></tr><tr><td></td><td></td><td>adjustments for color measurement</td></tr><tr><td></td><td></td><td>selection of address for IEEE488 interface</td></tr></table>	Menüpunkt	Untermenüpunkt	Funktion	Display Mode	Absolute	displays the measurement respective of any offset and calibration factors programmed	Relative Percent	displays the measurement as percentage of a reference value	Relative Factor	displays the measurement referred to a reference value	Relative Logarithmic Dose	displays the measurement in dB in relation to a reference value	No Display	accumulates the single readings and displays the result as exposure for measured quantity	Range		no display of the selected channel	Detector	Ampere	selects the range 1 to 8 or switches into the autorange mode (9)	Manual	shows the reading in ampere units	Select Table Entry	allows entering a calibration factor manually	Select Detector Recalibration	selects calibration data table entry	Offset		selects calibration data table entry	Reference		performs recalibration of selected calibration table entry with a reference lamp	Main Mode	CW	performs an automatic offset adjustment	Color	sets a reference value, used in different display modes	Pulse	measurement display of single channels	Pulse Color	measurement display color mode	RS232	measurement of short and single light pulses	Setup	IEEE488	measurement of short and single light pulses, result displayed as color values	Measurement Time	remote control RS232	Analog Out	remote control IEEE488	Dose Run Time	sets the integration time	Dose Maximum	selects the channel switched to analog output	Display Illumination	sets the maximum time duration for exposure measurement	Display Digits	sets the max. dose level for the exposure measurement mode	Color Setup	display backlight on / off		IEEE488 Address	adjusts the number of displayed digits (4, 5, or automatic)			adjustments for color measurement			selection of address for IEEE488 interface
Menüpunkt	Untermenüpunkt	Funktion																																																																				
Display Mode	Absolute	displays the measurement respective of any offset and calibration factors programmed																																																																				
	Relative Percent	displays the measurement as percentage of a reference value																																																																				
	Relative Factor	displays the measurement referred to a reference value																																																																				
	Relative Logarithmic Dose	displays the measurement in dB in relation to a reference value																																																																				
	No Display	accumulates the single readings and displays the result as exposure for measured quantity																																																																				
Range		no display of the selected channel																																																																				
Detector	Ampere	selects the range 1 to 8 or switches into the autorange mode (9)																																																																				
	Manual	shows the reading in ampere units																																																																				
	Select Table Entry	allows entering a calibration factor manually																																																																				
	Select Detector Recalibration	selects calibration data table entry																																																																				
Offset		selects calibration data table entry																																																																				
Reference		performs recalibration of selected calibration table entry with a reference lamp																																																																				
Main Mode	CW	performs an automatic offset adjustment																																																																				
	Color	sets a reference value, used in different display modes																																																																				
	Pulse	measurement display of single channels																																																																				
	Pulse Color	measurement display color mode																																																																				
	RS232	measurement of short and single light pulses																																																																				
Setup	IEEE488	measurement of short and single light pulses, result displayed as color values																																																																				
	Measurement Time	remote control RS232																																																																				
	Analog Out	remote control IEEE488																																																																				
	Dose Run Time	sets the integration time																																																																				
	Dose Maximum	selects the channel switched to analog output																																																																				
	Display Illumination	sets the maximum time duration for exposure measurement																																																																				
	Display Digits	sets the max. dose level for the exposure measurement mode																																																																				
	Color Setup	display backlight on / off																																																																				
	IEEE488 Address	adjusts the number of displayed digits (4, 5, or automatic)																																																																				
		adjustments for color measurement																																																																				
		selection of address for IEEE488 interface																																																																				
Sonstiges																																																																						
Bedieneinheit	16 Tasten, menügeführt																																																																					
Schnittstelle	RS232: 600 – 57600 Baud, 8 Data Bit, 1 Stop Bit, keine Parität, DSUB9 Stecker (weiblich) IEEE488: AH1, SH1, L4, T4																																																																					
Abmessungen	280 mm x 250 mm x 70 mm																																																																					
Garantie	12 Monate																																																																					

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Gigahertz Optik GmbH