XD-45-HUV

<https://www.gigahertz-optik.com/es-es/producto/xd-45-huv/>

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Descripción

The XD-45-HUV irradiance detector is designed for the evaluation of light exposure hazard values for artificial light sources. The three sensor design of this unique device covers the requirements for skin and eye risk assessment. It supports a number of standards and regulations:

- IEC 62471:2006 and EN 62471:2008
- 2006/25/EC
- EN 14255-1

ICNIRP/ACGIH spectral responsivity

The spectral effective function required for skin and eye risk evaluation is formed using two filtered sensors. This prevents the cross-talk and limited sensitivity between the UV-A, UV-B and UV-C spectral ranges inherent in a single filtered sensor solution. Using the two sensor method also produces a much better simulation of the intended spectral effectiveness function.

Additional radiometric UV-A responsivity

An additional UV-A sensor for the evaluation of UV-A_{315nm-400nm} human eye risk is included. All three sensors are mounted behind one 20 mm diameter cosine diffuser. For measurements of eye dependent radiance values a front adapter is supplied to limit the detector field-of-view to 80 degrees.

Traceable calibration

Calibration of the detector ICNIRP (W/m²) and UV-A (W/m²) responsivity is performed by the Gigahertz-Optik GmbH calibration laboratory for optical radiation measurements quantities. As with all light detectors supplied by Gigahertz-Optik calibration of absolute detector responsivity as well as detector individual measured relative spectral responsivity data is included.

Recommended optometer


For applications involving steady state light the X1₃ optometer is recommended. For unstable or short event light the optometer P-9801 should be selected.

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
General

Detector for UV Hazard



XD-45-HUV - UV-CB_{eff} and UV-A_{eff} Sensor - Typical Spectral Responsivity

XD-45-HUV - UV-A Sensor - Typical Spectral Responsivity

General	SUV(λ) IEC/EN 62471 spectrum detector useful for evaluation of artificial light sources. Set-up with in sum three photodiodes. Two are forming with a optical filter correction the SUV(λ) actinic UV-A und UV-BC spectral responsivity. A third detector with a radiometric UV-A responsivity is used for the measurement of the UV-A exposure of the human eye.
Especificación	
Campo de visión	20 mm diameter diffuser with removable 80 ° F.O.V. front adapter
Rango de medición	SUV(λ) (Σ of 2 cells): effective irradiance 0.5 mW/m ² to 10 W/m ² (max. resolution 0.05 mW/m ²), details see spectral plot.
Rango de medición	SUV (230 - 320) nm: effective irradiance 0.5 mW/m ² to 10000 W/m ² (max. resolution 0.05 μ W/m ²), details see spectral plot (usage starts at 230 nm).
Rango de medición	SUV (320 - 400) nm: effective irradiance 0.5 μ W/m ² to 10 W/m ² (max. resolution 0.05 μ W/m ²), details see spectral plot.
Rango de medición	UV-A (315 - 400) nm: irradiance 0.2 mW/m ² to 10000 W/m ² (max. resolution 0.02 mW/m ²), details see spectral plot.

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(remark: limited responsivity on the edge of the measurement range)

Miscelánea	
Peso	XD-45-HUV detector: 45 mm diameter, 30 mm high, 200 g (with cable) FOV adapter: 60 mm diameter, 32/49 mm high, 50 g
Nota	Specifications stated are for a X1-3 meter including the detector *) <i>The maximum measurable signal value may be limited by heat radiated by the source under test</i>
Garantía	12 month

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Gigahertz Optik GmbH (Sede central)
Gigahertz-Optik, Inc. (Oficina de EE.UU.)