


S-SDK-P9710

<https://www.gigahertz-optik.com/en-us/product/s-sdk-p9710/>

Product tags:




Description

Software Development Kit

Gigahertz-Optik offers programmers the S-SDK-P9710 Software Development Tool. The SDKs simplify integration of the P9710 and its variants in own software.

The S-SDK-P9710 is a C library (DLLs) that can be integrated into many programming environments (C#, LabView, Python, Matlab, .NET, Delphi, etc.). Beside the .dll's and LabView VI's a comprehensive documentation is delivered.

The programmer is able to control all functions of the device like: measurement routines, complex evaluations, etc.



Software Development Kit


Specifications

General

System requirements

- Minimum HDD space: 10MB (only files)
- Operating System: Windows 11, Windows 10 32-bit, Windows 10 64-bit, Windows 7 32-bit, Windows 7 64-bit
- communication: RS-232

Configurable with

Product Name	Product Image	Description	Go to product
P-9710		High-End Optometer for Measurement of CW-, Single Pulse and Modulated Radiation	https://www.gigahertz-optik.com/en-us/product/p-9710/

Purchasing information

Article-Nr	Modell	Description
Product		
15298231	S-SDK-P9710	Software Development Kit for the implementation of the P9710 or variants into custom made software

Contact, Calibration, Service & Support

We are known worldwide for excellent technical consulting and after sales support. Contact us to find together the best solution for you. Our services:

- Technical Consulting & Sales
- After-Sales Support
- Calibrations & Re-Calibrations ([ISO/IEC 17025 Calibration Services](#), [factory calibration](#), [Calibration of Third-Party Products](#))
- Repairs & Updates
- OEM & Feasibility Consulting of Customized Solutions

[Send us your inquiry](#) or contact us by phone or e-mail. We would welcome your feedback too or review us on [Google](#).

Gigahertz Optik GmbH (Headquarter)

Tel.: +49 (0)8193-93700-0
Fax: +49 (0)8193-93700-50
info@gigahertz-optik.de

An der Kaelberweide 12
82299 Tuerkenfeld, Germany

Gigahertz-Optik, Inc. (US office)

Phone: +1-978-462-1818
info-us@gigahertz-optik.com

Boston North Technology Park
Bldg B - Ste 205
Amesbury, MA 01913 USA