

CP-LG Series

<https://www.gigahertz-optik.com/zh-cn/product/cp-lg-series/>

产品标签:



Light Guide Accessory Series

The CP-LG Series offers all kind of optical light guides. We have many different length, detector designs (inline, 90°, probes, etc.) available.

Customized Solutions possible

For many applications customized designs are needed. This is starting by a different fiber length but also includes very sophisticated detector heads or Y-Fibers (2 in 1) or e.g. 1 to 3 detector heads if e.g. three BTS2048 shall measure with the same detector. We have a solution for this needs and provide many customized solutions as well.

Calibration of Light Guide

We also offer special calibrations for light guide solutions if needed. Our ISO 17025 calibration laboratory is very experienced in high quality calibrations.



BTS2048-xx-F with CP-CD-90-10 90° detector head



BTS2048-VL-F-Z01: Probe header for small LED chips (waferprober)

可配置的有

产品名称	产品图片	描述	转到产品
BTS2048-VL-F		Fiber-Coupled Version of the Versatile High Speed BTS2048-VL	https://www.gigahertz-optik.com/zh-cn/product/bts2048-vl-f/
BTS2048-VL-TEC-F		Fiber-Coupled Version of the Versatile High Speed BTS2048-VL-TEC	https://www.gigahertz-optik.com/zh-cn/product/bts2048-vl-tec-f/
BTS2048-UV-F		Measure in the UV with this high-quality compact fiber coupled spectroradiometer	https://www.gigahertz-optik.com/zh-cn/product/bts2048-uv-f/

联系、校准、服务和支持/Contact, Calibration, Service & Support

我们以出色的技术咨询和售后支持而闻名于世。请与我们联系，共同寻找适合您的最佳解决方案。我们的服务。

- 技术咨询和销售
- 售后支持
- 校准和重新校准（[ISO/IEC 17025校准服务](#)，工厂校准，[第三方产品的校准](#)）。
- 维修和更新
- 定制解决方案的OEM和可行性咨询

[请将您的询问发送给我们](#)

（*为必填项）或通过电话或电子邮件与我们联系。我们也欢迎你的反馈或在谷歌上评论我们。 [Google](#)

Gigahertz Optik GmbH (总部)

Tel.: +49 (0)8193-93700-0

Fax: +49 (0)8193-93700-50

info@gigahertz-optik.de

An der Kaelberweide 12

82299 Tuerkenfeld, Germany